Blood Glucose Test Strips

(Model: VGS01)

Package Insert

PRINCIPLE AND INTENDED USE

The **Keto-Mojo** Blood Glucose Test Strips (VGS01) are thin strips with a chemical reagent system. They work with the **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter as a system to quantitatively measure the glucose concentration in fresh capillary whole blood from fingertip. The **Keto-Mojo** blood glucose test is based on measurement of electrical current caused by the reaction of the glucose with the reagents on the electrode of the test strip. The blood sample is pulled into the tip of the test strip through capillary action. Glucose in the sample reacts with glucose enzyme and the mediator. Electrons are generated, producing a current that is positive correlation to the glucose concentration in the sample. After the reaction time, the glucose concentration in the sample is displayed. The meter is calibrated to display plasma-like concentration results.

The **Keto-Mojo** Blood Glucose Test Strips and **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter are intended for use outside the body (in vitro diagnostic use) by people with diabetes at home, as an aid to monitor the effectiveness of diabetes control. The system is for self-testing, intended to be used by a single person and should not be shared. The system should not be used for diagnosis of or screening for diabetes and not for use on neonates.

COMPOSITION

Each test strip (VGS01) contains the following reactive chemicals: Glucose oxidase < 25 IU, Mediator < 300 μ g.

Each test strip vial contains a drying agent. Each test strip pouch contains a drying agent.

STORAGE AND HANDLING

- Store test strips in a cool, dry place between 2-30°C (36-86°F). Keep away from heat and direct sunlight. Exposure to temperature and / or humidity outside the required condition may result in inaccurate readings.
- Do not freeze or refrigerate.
- Use the test strips at temperatures between 5-45°C (41-113°F).
- Use the test strips between 10-90% humidity.
- Do not store the meter, the test strips or control solution near bleach or cleaners that contain bleach.
- Replace the vial cap and close it tightly immediately after removing a test strip.
- Always keep the test strips in the original vial. Use the test strip immediately after removing it from the vial or foil pouch.
- Do not use your test strips beyond the expiry date (printed on the strip vial label or on the foil pouch or on the box) or discard date, whichever comes first, because this may cause incorrect test results.

Note: All expiration dates are printed in Year-Month format. 2024-01 indicates January, 2024.

- Write the discard date (6 months after first opening the vial) on the vial label when you first open it. Discard any remaining test strips after the discard date.
- Do not use test strips that are torn, bent, or damaged in any way. Do not reuse test strips.
- Keep the test strip away from children. Do not swallow test strips.
- Never ignore symptoms or make significant change to your diabetes control program without speaking to your healthcare professional.

PERFORMING A BLOOD GLUCOSE TEST

Materials provided: **Keto-Mojo** Blood Glucose Test Strips and package insert
Materials required but not provided: **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter, User's Manual, lancing device, clear cap (optional) and a new sterile lancet.
Refer to your User's Manual for complete instructions for blood sample collection before use.

- Wash your hands in with warm water and soap, dry them thoroughly.
- Prepare the lancing device.
- Check the expiration date (printed on the strip vial label or on the foil pouch or on the box). Do not use test strips beyond the expiration date.
- Insert the test strip into the meter. The meter turns on.
- Using a lancing device and new lancet to obtain a round drop of blood.
- Touch the blood drop to the strip tip until the meter beeps. Do not apply blood on the top of test strip.
- Your blood glucose test result will appear after the meter counts down from 5 to 1.

EXPECTED DIABETES CONTROL GOAL

Blood glucose values will vary depending on food intake, medication dosages, health, stress, or exercise.

The American Diabetes Association suggests the following targets for most non-pregnant adults with diabetes. More or less stringent glycaemic goals may be appropriate for each individual. Consult your health care professional for the target value that is appropriate for you.

Expected blood glucose levels for most non-pregnant adults with diabetes:¹

Time	Range, mmol/L
Before a meal	4.4 – 7.2
1-2 hours after beginning of the meal	Less than 10

Questionable or Inconsistent Results: Repeat the test according to the procedures described in your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter User's Manual. If attempts to correct a problem fail, contact your health care professional.

CHECKING THE SYSTEM

Use only **Keto-Mojo** Control Solutions (VGC01).

For complete details about checking the system, refer your User's Manual.

When to check:

- At least once a week
- When you open a new box of test strip
- When you want to check the meter and test strips
- If your test strips were stored in extreme temperature or humidity
- After cleaning your meter
- If you have dropped the meter
- Your test result does not match with how you feel

For confirmation of results, Control Solution Normal tests should fall within the **CTRL 2** range, and Control Solution High tests should fall within the **CTRL 3** range. When testing with Control Solution Normal, make sure you are matching the results to the **CTRL 2** range printed on the strip box or labels.

CAUTION: If your quality control test result falls outside the control range shown on the strip box or labels, **DO NOT** use the system to test your blood, as the system may not be working properly. If you cannot correct the problem, contact **Keto-Mojo** for further assistance.

LIMITATIONS

- The **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter, **Keto-Mojo** Blood Glucose Test Strips (VGS01) and Control Solution (VGC01) have been designed, tested and proven to work together effectively to provide accurate blood glucose measurements. Do not use components from other brands.
- Use only with whole blood. Do not use with serum or plasma samples.
- Only fingertip sampling is acceptable.
- Very high (above 70%) and very low (about 20%) hematocrit levels can cause false results. Talk to your health care professional to find out your hematocrit level.
- Abnormally high levels (above 3 mg/dL) of vitamin C and other reducing substances will produce false high blood glucose measurements.
- The system is tested to accurately read the measurement of glucose in whole blood within the range of 1.1 to 33.3 mmol/L.
- Fatty substances (triglycerides up to 1500 mg/dL or cholesterol up to 500 mg/dL) have no major effect on blood glucose test results.
- The **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System has been tested and shown to work properly up to 10,000 ft (3,048 meters).
- Severely ill persons should not run the glucose test with the **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.
- Patients using oxygen therapy are not recommended for testing with **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.
- Blood samples from patients in shock, or with severe dehydration or from patients in a hyperosmolar state (with or without ketosis) have not been tested and are not recommended for testing with **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.
- Dispose of blood samples and materials carefully. Treat all blood samples as if they are infectious materials. Follow proper precautions and obey all local regulations when disposing of materials.

PERFORMANCE CHARACTERISTICS

The **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter is calibrated by using YSI (Model 2300 STAT PLUS) Glucose Analyzer reference instrument, which is traceable to NIST reference standard.

Repeatability, Precision

Repeatability-Blood		
Interval	Glucose concentration	Standard Deviation (SD) or Coefficient of Variation (CV)
1	2.2 mmol/L	0.078 mmol/L
2	4.3 mmol/L	0.111 mmol/L
3	7.4 mmol/L	1.9%
4	10.8 mmol/L	2.2%
5	19.7 mmol/L	2.1%
Intermediate Precision-Control Solution		
Interval	Glucose concentration	Standard Deviation (SD) or Coefficient of Variation (CV)
1	2.2 mmol/L	0.078 mmol/L
2	6.7 mmol/L	2.1%
3	19.4 mmol/L	2.0%

System Accuracy

The capillary blood glucose measurements from 352 participants were taken by lay persons using the **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter with 3 strip lots (y). Capillary blood samples were obtained from fingertip for the **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter testing. Fingertip samples from the same subjects were also analyzed with YSI Model 2300 STAT PLUS Glucose Analyzer (x). The results were compared.

Linear Regression Results: Blood Glucose and β -Ketone Meter (y) vs. YSI Reference (x)				
Sample Site	Slope	Intercept (mmol/L)	R	N
Fingertip	0.9960	-0.0152	0.9974	352

Fingertip samples were used for YSI reference measurement.

The sample range was 2.1 to 30.49 mmol/L for **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter testing with blood sampled from fingertip sites.

Fingertip Site: System Accuracy Results for Glucose Concentration \geq 4.2 mmol/L			
Within \pm 5%	Within \pm 10%	Within \pm 15%	Within \pm 20%
184/285 (64.6%)	275/285 (96.5%)	285/285 (96.5%)	285/285 (96.5%)
Fingertip Site: System Accuracy Results for Glucose Concentration < 4.2 mmol/L			
Within \pm 0.28 mmol/L	Within \pm 0.56 mmol/L	Within \pm 0.83 mmol/L	
52/67 (77.6%)	67/67 (100%)	67/67 (100%)	

Fingertip Site: System Accuracy Results for Glucose Concentration			
Within \pm 5%	Within \pm 10%	Within \pm 15%	Within \pm 20%
217/352 (61.6%)	336/352 (95.5%)	352/352 (100%)	352/352 (100%)

For complete instructions, please refer to the User's Manual included with your meter. For additional questions or issues with this product, please contact **Keto-Mojo GK+** for further assistance.

REFERENCES

- ADA Clinical Practice Recommendations, 2021.

INDEX OF SYMBOLS

	Consult instructions for use		Use by		Contains sufficient for <n> tests
	For in vitro diagnostic use only		Lot number		Control range
	Temperature limitations		Manufacturer		Catalog number
	Use within 6 months of opening				Do not reuse

Keto-Mojo Canada LLC
700-401 West Georgia Street
Vancouver BC V6B 5A1

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