GK+ Blood Glucose and β -Ketone Dual Monitoring System

Dear **Keto-Mojo** System User,

Thank you for choosing the **Keto-Mojo GK+** system! The **Keto-Mojo GK+** System is designed for easy testing and control of blood Glucose and blood Ketone levels.

Read this User's Manual carefully before you use your meter system. Please keep your User's Manual in a safe place; you may want to refer it in the future. Please keep the meter box (and any strip boxes you purchase) as your control solution test ranges are printed or stickered on the box. You will also find manuals and other resources, such as instructional videos, at keto-mojo.com.

Thank you again for choosing the **Keto-Mojo GK+** System.

Principle and Intended Use

The **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System is designed to quantitatively measure the blood Glucose concentration and blood Ketone concentration in fresh capillary whole blood from fingertip. The **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System is based on measurement of electrical current caused by the reaction of the Glucose/Ketone with the reagents on the electrode of the test strip. The blood sample is pulled into the tip of the test strip through capillary action. Glucose/Ketone in the sample reacts with the enzyme and the mediator. Electrons are generated, producing a current that is positive correlation to the Glucose/Ketone concentration in the sample. After the reaction time, the Glucose/Ketone concentration in the sample is displayed. The meter is calibrated to display plasma calibrated results.

The **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System is intended for use outside the body (*in vitro* diagnostic use) by people with diabetes at home, as an aid to monitor the effectiveness of diabetes control. The system is for self-testing, intended to be used by a single person and should not be shared. The system should not be used for diagnosis of or screening for diabetes and not for use on neonates.

TABLE OF CONTENTS

1.	Understanding Your Testing Tools.....	1
	Your Meter System Overview.....	1
	Your Meter Display.....	2
	Important Safety Information.....	4
2.	Setting Up Your System.....	5
	Set the Date and Time.....	5
	Set the Beep Feature.....	7
	Set the Test Alarm.....	7
	Set the Meal Marker.....	8
	Set the Hypoglycemia (Hypo) Alarm.....	8
	Set the Ketone Warning.....	9
	Synching Your Meter with the Smartphone.....	10
3.	Performing a Test.....	11
	Preparing the Test Strip.....	11
	Preparing the Lancing Device.....	11
	Getting a Blood Drop and Testing.....	13
	Testing with Control Solution.....	16
	Using the Meter Memory.....	18
4.	Maintenance and Troubleshooting.....	19
	Replace the AAA Batteries.....	19
	Caring for Your Keto-Mojo GK+ Blood Glucose and β -Ketone Dual Monitoring System.....	20
	Troubleshooting Guide.....	22
5.	Technical Information.....	23
	System Specifications.....	23
	Limitations.....	24
	Warranty.....	24
	Index of Symbols.....	25

CHAPTER 1: UNDERSTANDING YOUR TESTING TOOLS

Your Meter System Overview

The Keto-Mojo GK+ Blood Glucose and β -Ketone Meter and Test Strips

Your Meter Display

The picture below shows all the symbols that appear on your meter display. With the meter off, press and hold the main button to see the complete display. To make sure the display is working properly, all of the segments should be clear and exactly like the picture below. If not, contact **Keto-Mojo** for further assistance.

Icon	What it Means
88/88	Top left area on the screen indicates date.
88:88	Top right area on the screen indicates year or time.
AM PM	The top middle area on the screen indicates morning or afternoon time.
d/m/d	Indicates the display form of date and month (d/m) or month and day (m/d).
Bluetooth symbol	When flashing indicates trying to connect. Solid icon is successful Bluetooth communication.
Bluetooth symbol with X	Indicates failed Bluetooth communication.
Alarm clock icon	Indicates alarms have been set.
GLU	Blood Glucose test mode.
KET	Blood Ketone test mode.
A	Indicates average value.
Battery icon with -	Indicates low battery or battery needs to be replaced.

MEM	Indicates test result history.
888	Center area on the display that shows test results or error codes.
	When strip is inserted the drop will flash, indicating the system is ready to test.
Ketone?	Ketone warning. The symbol of "ketone?" will appear when the blood glucose test result is equal to or above 16.7 mmol/L to indicate user perform blood ketone test to avoid ketoacidosis.
Hypo	Indicates that a low Glucose test result may cause hypoglycemia.
	Pre-meal marker.
	Post-meal marker.
	Control test result.
mmol/L	Test results are displayed as mmol/L according to local government regulation.
	Indicates the temperature is not suitable for testing.

Notes:

Your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter is pre-set with beep sound function, the meter will beep when:

- You press the main button.
- While setting the test type.
- When setting the date and time (in set-up mode).
- When the test strip is inserted and is ready for application of blood or control solution.
- When sufficient blood or control solution is pulled into the test strip.
- When the test is complete.
- If any error occurs during operation.

Meter Use and Precautions

- The meter is pre-set to display blood Glucose and blood Ketone concentration in millimoles per liter (mmol/L). This unit of measure cannot be adjusted.
- Do not get water or other liquids inside the meter.
- Keep the strip port area clean.
- Keep your meter dry and avoid exposing it to extremes in temperature or humidity. Do not leave it in your car.
- Do not drop the meter or get it wet. If you do drop the meter or get it wet, check the meter by running a quality control test. Refer to **Quality Control Test** for instructions.
- Do not take the meter apart. Taking the meter apart will void the warranty.

- Refer to the **Caring for Your Meter** section for details on cleaning the meter.
 - Keep the meter and all associated parts out of reach of children.
- Note:** Follow proper precautions and all local regulations when disposing of the meter and used batteries.

All systems preventive warnings with regard to EMC

1. This instrument is tested for immunity to electrostatic discharge as specified in IEC 61000-4-2. However, use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets, etc.) may cause damaging static discharges that may cause erroneous results.
2. This instrument complies with the emission and immunity requirements described in EN61326-1 and EN61326-2-6. Do not use this instrument in close proximity to sources of strong electromagnetic radiation. It may interfere with proper operation of the meter.
3. For professional use, the electromagnetic environment should be evaluated prior to operation of this device.

Important Safety Information

- Do not use the meter if it is dropped into water or water splashes on it.
- Test strips and lancets are for single use only.
- Do not drop blood on the flat surface of the test strip.
- Check the expiration dates and discard dates on your test strips foil pouch / kit box and control solution bottle label.
- Use only **Keto-Mojo** Blood Glucose/Blood Ketone Test Strip with your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter.
- Use only **Keto-Mojo** Blood Glucose/Blood Ketone Control Solution with your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter.

Potential Biohazard

Do not share your products with other persons. All products or objects, which come in contact with human blood, even after cleaning, should be handled as if capable of transmitting viral disease.

CHAPTER 2: SETTING UP YOUR SYSTEM

Set the Date and Time

Note: Before your first use of the meter for testing, it is important to adjust the meter settings to set the correct date and time, ensuring that results stored in the memory are captured correctly. Follow manual set-up only if not using the App. Use the **main button** to confirm a setting. If you push past a setting you want to set, keep going forward (you cannot go back), to start the set-up process again. To re-start manual set-up, press and hold the main button until the meter shuts off.

1. Enter the setting mode and set the clock

With the meter off, press and hold the main button for more than 3 seconds and release - the meter will beep, then all meter display segments will appear (then disappear) and 24h (or 12h) will flash indicating set-up mode.

Clock format: you have two options for time format - 12h or 24h. Press the buttons on the left side of the meter up or down to adjust and set the clock for 12 or 24-hour mode, then press the main button twice to save your choice. The meter will prompt you to set the date.

2. Set the date


The year position will now flash on the display. Press the side button to adjust the year until the year setting is correct, then press the main button to save your choice. The meter will now prompt you to set the date format.

The display form of date and month will now flash, press the side button to set the date display form for m/d or d/m mode, then press the main button once to save your choice. The meter will prompt you to set the month.

The month will now flash, press the side button to adjust the month. Press the main button until the meter beeps to set.


The day will now flash, press the side button to adjust the date until the date setting is completed, then press the main button until the meter beeps to set.

3. Set the time

After the date setting is completed, press the main button until the meter beeps to set, the hour will now flash. Press the side button to adjust the current hour until the hour setting is completed, then press the main button until the meter beeps to set.

The minute will now flash, press the side button to adjust the minute until the minute setting is completed. Press the main button once to save your choice.

Note:

Before your first use of the meter for testing, please adjust the meter settings to set the correct date and time, ensuring that results stored in the memory are shown with the correct date and time.

Set the Beep Feature

After setting the time, press the main button until the meter beeps to set, the word "beep" will now flash on the screen. Press the side button to set ON or OFF, press the main button to save the beep setting.

**Set the Test Alarm**

After the audio setting is completed, the meter will enter the test alarm setting mode; this reminder function reminds users to take a measurement. You can set up to 5 reminders per day. If you turn 5 test alarms on (A1, A2, A3, A4 and A5), your meter is pre-set with the following times for your convenience, you can adjust each time to suit your needs.

A1 7:00 A2 9:00 A3 14:00 A4 18:00 A5 22:00

Before you set the time, the ⌚ and the word "OFF" will be displayed, symbol "A1" flashes. Press side button to turn the alarm function on or off, and press the main button to save.


If you select "On", the hour flashes, "A1" and the ⌚ will remain on the display. Press side button to select the hour. Press the main button once to set.

When the minute flashes, press the side button to select 00, 15, 30, or 45. These are the only choices. Press the main button once to set.

The next alarm "A2" flashes on the display with "OFF".

You can set a second alarm by pressing the side button to turn on the second alarm. Perform the same procedures to set the rest alarms.

Note: If the meter is on at the test reminder time, the test alarm will not be activated.

Set the Meal Marker

After setting the test alarm functions, the symbols of 🍎🍷 will now flash, along with the word "On" on the display. Press the side button to turn the meal marker function on or off, press the main button once to set

Set the Hypoglycemia (Hypo) Warning

After setting the meal marker, you can set the hypo alarm, which indicates a possible hypoglycemic condition (blood sugar level too low).

After you have confirmed the selection of setting test alarms, the **Hypo** flashes on the display along with "On" on the display. Press the side button to turn the Hypo alarm function on or off, then press the main button once to set. If you select the hypo alarm to be "on", the display shows 3.9 mmol/L with the symbol of "GLU" on the display, press the main button to set.

Notes:

- Talk to your healthcare professional to help you decide the hypo level that suits your health condition.
- For healthcare professionals, the hypoglycemic level may vary from person to person. It is recommended that you turn the hypo alarm function **OFF** when using the meter in a professional facility.


Set the Ketone Warning

After the hypo warning setting is completed, the **Ketone?** symbol will now flash, along with word “On” on the display. Press the side button to turn the Ketone Warning on or off; press the main button to set.

If you turn the Ketone Warning on, the **Ketone?** symbol will appear with the other symbols that you have set before.

After setting the Ketone Warning, the screen will show all symbols which you have set before. Press the main button and the meter will turn off.

Once all the settings are completed, if you want to change the setting, please press and hold the power main button for 3 seconds when the meter is powered off and then to return to the set-up mode.

Synching Your Meter with the Smartphone

Synching prepares your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter and Smartphone to communicate with each other. The distance between the meter and Smartphone should be within 5 meters. Download the **Keto-Mojo App** before synching your meter and Smartphone.

Note: The **Keto-Mojo App** is compatible with Android software version 8 and above and iOS version 11 and above.

WARNING

Do Not synch another person's meter with your Smartphone.

Synching allows your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter and Smartphone to communicate with each other. This sets your meter date and time to your Smartphone date and time (manual set-up will then not be necessary).

Follow the guidelines and instructions to successfully synch the meter with the **Keto-Mojo App**.

Please note that data and reports in the **Keto-Mojo GK** App are not used for medical care reference by user or healthcare professionals.

Instructions for pairing your meter with the **Keto-Mojo** iOS and Android App can be found by either scanning the QR code below, or by typing the following weblink into a new browser tab:

<https://keto-mojo.com/bluetooth-app-connection/>

CHAPTER 3: PERFORMING A TEST

Set up your meter correctly and have all the materials you will need ready before you begin testing. This includes your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter, **Keto-Mojo** Blood Glucose Test Strips, **Keto-Mojo** Blood Ketone Test Strips, and **Keto-Mojo** Lancing Device and Lancets.

Preparing the Test Strip

1. Wash and dry your hands well before testing.
2. Remove a test strip from the test strip foil pouch.
3. With the meter off, insert the test strip into the meter in the direction of the arrows. (there is no need to press the main button to begin testing).

4. A symbol with a test strip with a flashing blood will appear letting you know the meter is ready to test.


Note:

Check the expiration and discard dates on the test strip foil pouch. All expiration dates are printed in Year-Month format. 2024-01 indicates January, 2024. Make sure the test strip does not appear damaged. Prior to testing, wipe and dry the test site with an alcohol swab or soapy water. Make sure there is no cream or lotion on the test site.

Preparing the Lancing Device


For fingertip sampling, adjust the depth penetration to reduce the discomfort. You do not need the clear cap for fingertip sampling.

1. Unscrew the lancing device cover from the body of the lancing device. Insert a sterile lancet into the lancing device and push it until the lancet comes to a complete stop in the lancing device.

Note: The **Keto-Mojo** lancing device uses **ONLY Keto-Mojo** sterile lancets.

2. Hold the lancet firmly in the lancing device and twist the safety tab of the lancet until it loosens, then pull the safety tab off the lancet. Save the safety tab for disposing the used lancet.

3. Carefully screw the cover back onto the lancing device. Avoid contact with the exposed needle. Make sure the cover is fully sealed on the lancing device.

4. Adjust the puncture depth by rotating the lancing device cover. There are several different puncture depth settings. To reduce discomfort, use the lowest setting that still produces an adequate drop of blood.

Adjustment:

- 1 for delicate skin
- 2 and 3 for normal skin
- 4 and 5 for calloused or thick skin


Note: Greater pressure of the lancing device against the puncture site will also increase the puncture depth.

Getting a Blood Drop for Testing


1. Pull the cocking barrel back to set the lancing device. You may hear a click to indicate the lancing device is now loaded and ready for obtaining a drop of blood.

2. Press the lancing device against the side of the finger to be lanced with the cover resting on the finger. Push the release button to prick your fingertip. You should hear a click as the lancing device activates.

3. Wipe away the first drop of blood with a clean paper towel to ensure a more accurate result. Gently massage from the base of the finger to the tip of the finger to obtain the required blood volume (half the size of a match head). Avoid smearing the drop of blood. For the greatest reduction in pain, lance on the side of the fingertips. Test immediately after a good blood drop has formed.

4. Immediately touch the tip of the test strip to the drop of blood. The blood will be pulled into the test strip through the tip. Make sure that the blood sample has fully filled the check window of the tip of the strip. Hold the tip of the test strip in the blood drop until the meter beeps.


Note:

If the blood sample does not fill the check window, do not add a second drop. Discard the test strip and start over with a new test strip.

5. For Glucose testing, the meter counts down from 5 to 1, and for Ketone testing, the meter counts down from 9 to 1. The result then appears on the display. The test result will automatically be stored in the meter memory. Please do not touch the test strip during the countdown as this may result in an error.

Discard the Used Test Strip

You can eject and discard the used test strip by using the strip ejector. Meter turns off automatically after a beep.

Potential Biohazard

Dispose of the used test strips as medical waste.

Questionable or Inconsistent Results:

If your blood Glucose/blood Ketone result does not match how you feel, please:

- Check the expiration date of the test strip.
- Confirm that the room temperature during testing blood Glucose is between 5 and 45°C, and between 7.5 and 45°C for testing blood Ketone.
- Make sure the test strip has been stored in cool, dry place.
- Make sure the test strip was used immediately after removing from the test strip foil pouch or vial.
- Make sure that you followed the test procedure correctly.
- Perform a control solution test (See **Performing a Control Test** for instructions).

After checking all of the conditions listed above, repeat the test with a new test strip. If you are still unsure of the problem, please contact **Keto-Mojo** immediately for further assistance.

Removing the Used Lancet

Unscrew the lancing device cover. Place the safety tab of the lancet on a hard surface and carefully insert the lancet needle into the safety tab.

Press the release button to make sure that the lancet is in the extended position. Slide the ejection button forward to discard the used lancet. Place the lancing device cover back on the lancing device.

Potential Biohazard

Always dispose of the used lancet properly to prevent injury or contamination to others.

Caution:

- *Do not use the lancet if the safety tab is missing or loose when you take the lancet out of the bag.*
- *Do not use the lancet if the needle is bent.*
- *Be cautious whenever the lancet needle is exposed.*
- *Never share lancets or the lancing device with other people to prevent possible infections.*
- *In order to reduce the risk of infection from prior use of the instrument, always use a new, sterile lancet. Do not reuse lancets.*
- *Avoid getting the lancing device or lancets dirty with hand lotion, oils, dirt or debris.*

Testing with Control Solution

Why Perform Control Tests

Performing a control test lets you know that your meter and test strips are working properly to give reliable test results. Suggestions for when to perform a control test:

- At least once a week(optional)
- When you open a new box of test strip
- When you want to check the meter and test strips
- If your test strips were stored in extreme temperature or humidity
- After cleaning your meter
- If you have dropped the meter
- Your test result does not match with how you feel

About the Control Solutions

- Only use **Keto-Mojo** Glucose/Ketone Control Solutions (Level 2 and 3) to practice on the system.
- The control solution results are not included in the average value calculation.
- All expiration dates are printed in Year-Month format. 2022-01 indicates January, 2022.
- Do not use control solution that is past the expiry date or discard date (the control solution will expire 6 months after the bottle is opened for the first time).
- Shake the bottle well before use.
- Close the bottle tightly after use.

Performing a Control Test

1. Remove a test strip from the test strip foil pouch.

Note: Check the expiration of the test strips. Do not use expired test strip.


2. With the meter off, insert a test strip into the meter in the direction of the arrows. It is not required to press the main button to begin testing.

3. The meter turns on after a beep. An image of a test strip with a flashing blood drop will appear letting you know the meter is ready to test.

4. Shake the control solution bottle thoroughly. Squeeze the control solution bottle gently and discard the first drop. Squeeze out a second small drop on a clean nonabsorbent surface.

Note: Do not apply control solution to the test strip directly from the bottle.

5. Immediately touch the tip of the test strip to the drop of control solution. The control solution is pulled into the test strip through the strip tip.

Note: If the control solution sample does not fill the check window, do not add a second drop. Discard the test strip and start over with a new test strip.

6. Hold it in the drop until the meter beeps, and then you see the meter count down on the screen, followed by your control test result.

Note: The meter will automatically recognize and mark the control result for you.
Control results are not included in the 7, 14 and 30 day average calculation.

Understand Your Control Test Result

Compare your control test result with the ranges printed on the test strip box or label.

Notes: If your control test result is out of range:

- Check the expiration dates and discard dates of the test strip and control solution. Make sure that the control solution bottle has not been opened for more than 6 months. Discard any expired test strips or control solution.
- Confirm the temperature in which you are testing of Glucose control is between 10 and 40°C, and testing of Ketone control is between 15 and 40°C
- Make sure that the control solution bottle has been tightly capped.
- Make sure the test strip was used immediately after removing from the test strip foil pouch.
- Make sure the control solution was mixed well.
- Confirm that you are using **Keto-Mojo** control solution.

- *Make sure that you followed the test procedure correctly.*

After checking all of the conditions listed above, repeat the control solution test with a new test strip. If your results still fall out of the range indicated on the test strip box or label, your meter or test strips may not be working properly. **DO NOT** use the system to test blood. Contact **Keto-Mojo** for further assistance.

To turn your meter off, just remove the test strip. Dispose of the used test strips as medical waste. The result will be automatically marked and stored in the meter memory. Control results will not be included in your blood Glucose averages.

Using the Meter Memory

Your meter automatically stores up to 1000 results with time and date. Test results are stored from the newest to the oldest. The meter will also calculate the average values of blood Glucose records from the last 7, 14 and 30 days.

Notes:

- *If there are already 1000 records in memory, the oldest record will be erased to make room for a new one.*
- *It is very important to set the correct time and date in the meter, please make sure the time and date are correct after you change your battery.*
- *Control results of blood Glucose are not included in the 7, 14 and 30 day average calculation.*

Viewing Your Test Results

When your meter is off, press the main button to turn meter on. After a beep, a symbol of strip flashes on the display, press the main button again, the 7-day average of blood Glucose will appear in the center of the display. If you want to review the memory after you immediately performed a test, with the test result on the display, press the main button to see the 7-day of blood Glucose average.

Continue to press the main button to view the 14-day average of blood Glucose, then press the main button again to review the 30-day average of blood Glucose.

Continue to press the main button to review previous results in order (including blood Glucose and blood Ketone). Results will be shown starting with the most recent; each result will show the date and time the test was taken.

When **END** appears on the display, you have viewed all of the results in the memory.


CHAPTER 4: MAINTENANCE AND TROUBLESHOOTING

Proper maintenance is recommended for best results.

Replace the AAA Batteries

When the meter AAA batteries need to be replaced, the battery symbol () will appear, or you may receive an **E-11** error code.

1. Turn your meter off before changing the batteries.
2. Press firmly on the battery cover and slide.

Note:

After you change the AAA batteries, your meter prompts you to confirm the meter's time and date settings. All of the test results are saved in the memory.

3. Lift out and remove the old AAA batteries.
4. Place the new AAA batteries under the prongs and into the battery compartment.

5. Slide the battery cover back into place, lining up with the open slots, and close firmly.

Warning:

Keep batteries away from children. A lithium battery is poisonous. If swallowed, immediately contact your doctor or poison control center. Discard batteries according to your local environmental regulations.

Caring for Your Keto-Mojo GK+ Blood Glucose and β -Ketone Dual Monitoring System

- Store meter in the carrying case provided whenever possible.
- Wash and dry hands well after handling to keep the meter and test strips free of water and other contaminants.
- The Blood Glucose and β -Ketone Meter is a precision electronic instrument. Please handle it with care.
- Avoid exposing meter and test strips to excessive humidity, heat, cold, dust, or dirt.

Cleaning and Disinfection

Your **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter should be cleaned and disinfected a minimum of once per week. Use only Clorox™ Healthcare Bleach Germicidal Wipes, which has been proven to be safe to use with the **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter.

Cleaning is part of your normal care and maintenance and should be performed prior to disinfection, but cleaning does not kill germs. After use and exposure to blood, all parts of this kit can potentially transmit infectious diseases. Disinfecting reduces the risk of transmitting infectious diseases.

Note: If the meter is being operated by a second person who is providing testing assistance to you, the meter should be cleaned and disinfected prior to use by the second person.

1. Cleaning Your Meter

Step 1: Take one piece of Clorox™ Healthcare Bleach Germicidal Wipes (EPA Registration No. 67619-12) from the container.

Step 2: Clean the entire meter surface including front side, back side, right side and left side.

The meter should be cleaned when it is visibly dirty or at a minimum of once per week. This cleaning is to prepare the meter surface for a disinfection process.

2. Disinfecting Your Meter

Step 1: After cleaning your meter, take out a new piece of Clorox™ Healthcare Bleach Germicidal Wipes.

Step 2: Wipe the entire surface including front side, back side, right side and left side of the meter, with a back and forth movement.

Step 3: Keep the meter surface wet for at least one minute.

Step 4: Wait for the surface of meter to dry.

Clorox™ Healthcare Bleach Germicidal Wipes contain Sodium hypochlorite 0.55%, which has been proven to be safe to use with the Blood Glucose and β -Ketone Dual Monitoring System. Clorox™ Healthcare Germicidal Bleach Wipes are available by contacting Krasity's Medical Supply at 800-537-1394 directly or purchasing at <http://www.walmart.com> and <http://www.staples.com/>.

The meter should be disinfected a minimum of once per week. The meter disinfection process has been validated for 608 disinfection cycles, which is equivalent to cleaning and disinfecting your meter every 3 days for 5 years. This is to ensure that your meter will operate properly over the 5-year life of the meter.

Note:

- Do not use alcohol or any other solvent.
- Do not allow liquid, dirt, dust, blood, or control solution to enter the test strip port or the data port.
- Do not squeeze or wipe gauze into the test strip port.
- Do not spray cleaning solution on the meter.
- Do not immerse the meter in any liquid.

Note: Although it has not been observed, some alterations may appear on your meter due to the cleaning and disinfection procedure. Items such as: cloudy display window, plastic housing cracking, meter buttons that do not function, partial display on full screen, unable to execute the meter's initial set up, etc. may appear. If you notice any of these external changes to your meter or any changes to the performance of your meter stop using the meter and please contact the local distributor for further assistance. Please contact your healthcare professional if you need help.

If you have questions about cleaning or disinfection, or if you see evidence of physical damage, please contact the local distributor for further assistance. Please contact your healthcare professional if you need.

Troubleshooting Guide

What You See	What It Means	What You Should Do
E 1	Add sample or control solution before the flashing blood drop appears.	Discard the test strip and repeat the test with a new test strip. Add sample or control solution after you see the flashing blood drop on the display.
E 2	The meter is sensing the use of a used or contaminated test strip.	Discard the test strip and repeat the test with a new test strip. Wait until you see the flashing blood drop on the display before testing.
E 3	Incorrect test strip.	Discard the test strip and repeat the test with a new test strip. Make sure that you are using Keto-Mojo Glucose / Ketone test strip.
E 4	Incorrect sample.	Discard the test strip and repeat the test with a new test strip. Make sure that you are using a blood sample or Keto-Mojo Glucose/Ketone control solution.
E 5	Temperature out of range.	Move to an area that is within the operating temperature for the meter. Let the meter adjust to this temperature for 20 minutes before performing a test.
E 6 E 7	Potential hardware issue.	Take out batteries and restart the meter. If the problem continues, contact Keto-Mojo immediately for further assistance.
E 10	Insufficient sample.	Repeat test and apply enough sample to fill check window of the test strip.
E 11	Low Battery Indicator	Replace batteries and check date/time/year settings.
HI	Blood Glucose test result is above 33.3 mmol/L; blood Ketone test result is above 8.0 mmol/L.	Wash and dry your hands, then repeat the test using a new test strip. If your result still flashes HI , contact your health care professional as soon as possible.
LO	Blood Glucose test result is below 1.1 mmol/L; blood Ketone test result is below 0.1 mmol/L.	Wash and dry your hands, then repeat the test using a new test strip. If your result still flashes LO , contact your health care professional as soon as possible.

CHAPTER 5: TECHNICAL INFORMATION

System Specifications:

Feature	Specification
Measurement range	Blood Glucose: 1.1-33.3 mmol/L; Blood Ketone: 0.1-8.0 mmol/L
Result calibration	Plasma-equivalent
Sample	fresh capillary whole blood
Sample volume:	Blood Glucose: about 0.5 μ L; Blood Ketone: about 0.8 μ L
Test time	Blood Glucose: about 5 seconds; Blood Ketone: about 10 seconds;
Power source	Two AAA LR03 1.5V batteries
Battery life	12 months or approximately 1,000 tests
Units of measure	Blood Glucose: the meter is pre-set to either millimoles per liter (mmol/L); Blood Ketone: millimoles per liter (mmol/L) is the default setting
Memory	Up to 1000 records with date and time
Automatic shutoff	2 minutes after last action
Dimensions	90.4 mm x 54.5 mm x 27.8 mm
Display size	39 mm x 41 mm
Weight	Approximately 63g
Operating temperature	Blood Glucose: 5-45°C Blood Ketone: 7.5-45°C
Operating relative humidity	10-90% (non-condensing)
Hematocrit range	Blood Glucose: 20-70%; Blood Ketone: 20-65%;
Bluetooth	Bluetooth Low Energy (BLE)

Limitations

The **Keto-Mojo GK+** Blood Glucose and β -Ketone Meter, **Keto-Mojo** test strips and control solutions have been designed, tested and proven to work together effectively to provide accurate blood Glucose/blood Ketone measurements. Do not use components from other brands.

- Use only with whole blood. Do not use serum or plasma samples.
- Only fingertip sampling is acceptable.
- Hematocrit levels beyond the hematocrit range of the system can cause false results. Talk to your healthcare professional to find out your hematocrit level.
- Abnormally high levels of vitamin C and other reducing substances will produce false high blood Glucose measurements.
- The **Keto-Mojo GK+** system is tested to accurately read the measurement of Glucose in whole blood within the range of 1.1-33.3 mmol/L; Ketone in whole blood within the range of 0.1-8.0 mmol/L.
- The **Keto-Mojo GK+** Blood Glucose and β -Ketone dual Monitoring System has been tested and shown to work properly up to 10,000ft (3,048 meters) for blood Glucose testing and 8,700 ft (2,651 meters) for blood Ketone testing.
- Severely ill persons should not run the Glucose and Ketone test with the **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.
- Patients using oxygen therapy are not recommended for testing with **Keto-Mojo** Blood glucose test strips (VGS01).
- Blood samples from patients in shock, with severe dehydration or in a hyperosmolar state are not recommended for testing with **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.
- Dispose of blood samples and materials carefully. Treat all blood samples as if they are infectious materials. Follow proper precautions and obey all local regulations when disposing of materials.

Warranty

Please register your warranty online at keto-mojo.com.

Index of Symbols

These symbols may appear on the packaging and in the instructions for the **Keto-Mojo GK+** Blood Glucose and β -Ketone Dual Monitoring System.

Symbols	Explanation
	Consult instructions for use
	For <i>in vitro</i> diagnostic use only
	Manufacturer
	Temperature limitations
	Contains sufficient for <n> tests
	Use by
	Lot Number
	Sterilized using irradiation
	Control range
	Catalog number
	Model number
	Do not re-use
	Caution, consult accompanying documents
	Dispose items according to local relevant laws regarding disposal and recycle
	Use within 6 months of opening
	Serial Number
	Keep away from sunlight and heat

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