



ICON

DESCRIPTION

08/08

Top left area on the screen indicates date

AM 08:08  
PM 08:08

Top right area on the screen indicates the year or time

d/m/d

Date, month, date



Successful Bluetooth communication



Failed Bluetooth communication

GLU

Blood glucose test mode

KET

Blood ketone test mode

A

Indicates average value



Indicates low battery or they need to be replaced (2 AAA)

MEM

Indicates test result history

00.0

Center area on the display that shows the test results or error codes



Indicates the system is ready to test

Ketone?

Ketone warning

Hypo

Indicates that a low test result may cause hypoglycemia



Pre-meal marker



Post-meal marker

mg/dL

Blood glucose test results are displayed as mg/dL

mmol/L

Blood  $\beta$ -ketone test results are displayed as mmol/L



Indicates the temperature is not suitable for testing